




**CERVICAL DILATATION/
EFFACEMENT SIMULATORS
LF01069U
INSTRUCTION MANUAL**



Life/form® Products by NASCO

About the Simulators...

The **Life/form**[®] Cervical Dilatation/Effacement Simulators are designed to facilitate the teaching, demonstration, and practice of cervical examination prior to birth. Every effort has been made to provide exceptional realism in these simulators so that the visual and tactile sensations experienced provide a positive transition to patient procedures.

Internal Structure

The complete set of simulators consists of six separate pelvic blocks. The external appearance of all six blocks is identical with the variations being in the internal structure. Internal texture, tissue density, and correct anatomical size are all carefully represented to provide the most realistic condition possible.

The simulators depict six different cervical conditions prior to birth. Each pelvic block is clearly identified on the bottom to indicate which condition exists internally. The different conditions depicted are:

Stage 1	No dilation No effacement
Stage 2	2 cm dilation 50% effacement
Stage 3	2 cm dilation Full effacement
Stage 4	5 cm dilation Full effacement
Stage 5	7 cm dilation Full effacement
Stage 6	9 cm dilation Full effacement

With predetermined measurements of the internal cervical condition, all guesswork or estimating on the part of the student and instructor is eliminated. Everyone's tactile understanding of 2, 5, 7, and 9 cm dilation and effacement condi-

tions is identical, eliminating misconceptions and inaccurate patient condition evaluations.

General Instructions for Use

The use of gloves and liberal use of lubricant teach and reinforce good practices and are essential to protect the physical properties of the simulators. The simulators must be treated like patients. Rough usage and lack of lubrication will generate tears in the soft tissue walls.

Care of the Simulators

Store simulators in clean plastic bags. The soft plastic used in the simulators does attract dust. As in a real situation, the user should practice care and good hygiene to assure the lasting life and clean appearance of the simulators.

Solvents or other corrosive materials will damage the simulators. Never place simulators on newsprint or any kind of printed paper or plastic. Avoid contact with ball point pens. These materials will transfer an indelible stain.

Normal soil accumulated on the surface of the simulators can be removed with mild soap and lukewarm water. Use REN Cleaner (W9919UA) to remove stubborn stains from simulators. Simply spray the soiled area and wipe clean with a soft cloth or paper towels.

Cautions

DO NOT use any petroleum based lubricant. Such materials will cause rapid deterioration and permanent damage. Use Lubricant Kit (LF00985US) or Ivory liquid detergent only.

Supplies/Replacement Parts

W09919UA	REN Cleaner
LF00985UA	Lubricant Kit

Other Available *Life/form* Simulators

- LF00698U** Adult Injectable Arm (White)
LF00856U Female Catheterization
LF00901U Prostate Examination
LF00906U Ostomy Care
LF00929U Surgical Bandaging
LF00957U Enema Administration
LF00958U Pediatric Injectable Arm
LF00961U Intramuscular Injection
LF00984U Breast Examination
LF00995U Arterial Puncture Arm
LF00997U Adult Injectable Arm (Black)
LF00999U Pediatric Injectable Head
LF01008U Intradermal Injection Arm
LF01012U Heart Catheterization (TPN)
LF01019U Ear Examination
LF01020U Supplementary Ear Set
LF01025U Male Cath-Ed I
LF01026U Female Cath-Ed II
LF01027U Peritoneal Dialysis
LF01028U Suture Practice Arm
LF01036U Spinal Injection
LF01037U Hemodialysis Arm
LF01062U Pelvic, Normal & Abnormal
LF01063U Stump Bandaging, Upper
LF01064U Stump Bandaging, Lower
LF01069U Cervical Effacement
LF01070U Birthing Station
LF01082U Cricothyrotomy
LF01083U Tracheostomy Care
LF01084U Sigmoidoscopic Examination
- LF01087U** Central Venous Cannulation
LF01095U Blood Pressure Arm
LF01108U Intraosseous Infusion Simulator
- LF01121U** Advanced IV Arm
LF01139U Advanced IV Hand
LF01142U Auscultation Trainer
LF01162U Venatech IV Trainer
LF03000U **CPARLENE®** Series
LF03601U Adult Airway Management Trainer
LF03602U Adult Airway Management on Manikin
LF03603U Adult Airway Management Head Only
LF03609U Child Airway Management Trainer
- LF03610U** Child Airway Management Trainer Head Only
LF03611U Child Defibrillation Chest Skin
LF03612U Child IV Arm
LF03613U Child Blood Pressure Arm
LF03614U Child Intraosseous Infusion/Femoral Access Leg Only
LF03615U Complete Child **CRiSis™** Update Kit
LF03616U Child **CRiSis™** Manikin
LF03617U Deluxe Child **CRiSis™** Manikin with Arrhythmia Tutor
LF03620U PALS Update Kit
LF03621U Infant Airway Management Trainer Head Only
LF03622U Intraosseous Infusion Right Leg
LF03623U Infant Airway Management Trainer
LF03626U Child Femoral Access Injection Pad Replacement
LF03632U Child Intraosseous Infusion/Femoral Access Leg on a Stand
LF03633U Child Airway Management Trainer with Torso
LF03693U **Basic Buddy** CPR Manikin
LF03699U "Airway Larry" Airway Management Trainer
LF03720U **Baby Buddy** Infant CPR Manikin
LF03953U **CRiSis™** Manikin
LF03955U Deluxe **CRiSis™** Manikin
LF03965U Deluxe "Plus" **CRiSis™** Manikin
LF04001U **GERi™** Nursing Manikin
LF04020U **KERi™** Nursing Manikin
LF04021U **KERi™** Basic Manikin
LF04022U **KERi™** Advanced Manikin
LF04030U **GERi™** Advanced Manikin
LF04040U **GERi™** Basic Manikin
LF06001U **CPR Prompt™** Adult/Child Manikin
LF06012U **CPR Prompt™** Infant Manikin
LF06200U **CPR Prompt™** Keychain Rescue Aid
LF06200U **CPR Prompt™** Rescue and Practice Aid

Nasco Fort Atkinson

901 Janesville Avenue, P.O. Box 901
Fort Atkinson, Wisconsin 53538-0901
1-800-558-9595

www.eNASCO.com • E-mail: lifeform@eNASCO.com